


MT9520 Binocular**MT9530 Trinocular****MT9540 Ergonomic Binocular**

- SWH10X Widefield high eyepoint eyepiece, FN22, tube OD 30mm
- Strain Free PlanAchromat Objective Lenses, Infinity Corrected
 - 4X, NA 0.10, WD 25.89mm
 - 10X, NA 0.25, WD 7.3mm
 - 40X, NA 0.65, WD 0.72mm
 - 20X, NA 0.40, WD 5.0mm (optional)
 - 60X, NA 0.85, WD 0.2mm (optional)
 - 100X oil, NA 1.25, WD 0.14mm (optional)
- Abbe Strain Free Condenser, NA 1.25, with iris diaphragm
- Ceramic coated flat-top stage, 190mm x 140mm, with 80mm x 52mm X-Y movement; low-positioned, drop down controls
- Slide specimen holder
- Intermediate tube with sliding analyzer
- Rotatable polarizer with full wave plate, attached below condenser
- Daylight balancing filter
- Dust cover

