


## MT9600 Series Phase Contrast Gout Microscopes

### MT9620 Binocular

### MT9630 Trinocular

### MT9640 Ergonomic Binocular



- SWH10X Widefield high eyepoint eyepiece, FN22, tube OD 30mm
- U Plan Phase Achromat Objective Lenses, Infinity Corrected
  - 10X, NA 0.25, WD 10.7mm
  - 20X, NA 0.40, WD 7.29mm
  - 40X, NA 0.65, WD 0.5mm
  - 100X oil, NA 1.25, WD 0.23mm
- Zernicke Phase Turret Condenser with stops for Ph1, Ph2, Ph3, Darkfield (up to 40X lens), Brightfield with iris diaphragm
- Ceramic coated flat-top stage, 190mm x 140mm, with 80mm x 52mm X-Y movement; low-positioned, drop down controls
- Slide specimen holder
- Intermediate tube with sliding analyzer
- Rotatable polarizer with full wave plate, attached below condenser
- Daylight balancing filter
- Dust cover



Phase Turret Condenser



Polarizer in Swing out mount

